

Urine Control With Biochemical Diagnostics

Detectabuse™ AU/NZ MIX For Onsite Devices



Top Reasons to use Detect abuse™ Controls:

- Human urine-based matrix simulates patient samples
- Recommended by the majority of handheld device manufacturers
- 12 drug constituents within one control vial
- Stable for up to 3 years frozen, or refrigerated (unopened)
- Ready to use in 5 mL vial sizes
- Use Detect abuse™ controls to train staff in the accurate use of hand-held devices
- Use Detect abuse™ controls to evaluate the performance of various devices on the market
- Each Detect abuse™ control contains: THC, Benzoylcegonine (COC), PCP, Morphine (OPI), Methamphetamine, Amphetamine, Secobarbital (BAR), Oxazepam (BZO), Methadone, Methaqualone, Propoxyphene, Nortriptyline (TCA).

The Detect abuse Advantage:

Detect abuse AU/NZ MIX controls are formulated to work with all known handheld drug-screening devices.

NEW! AUSTRALIAN-NEW ZEALAND CUTOFF LEVELS

TARGET VALUES	-50%	Cutoff*	+50%
Delta-9-THC-COOH	25	50	75
Benzoylcegonine	150	300	450
Phencyclidine (PCP)	12.5	25	37.5
Morphine	150	300	450
d-Methamphetamine	150	300	450
d-Amphetamine	150	300	450
Secobarbital	150	300	450
Oxazepam	100	200	300
Methadone	150	300	450
Methaqualone	150	300	450
Propoxyphene	150	300	450
Nortriptyline	500	1000	1500

Detectabuse AU/NZ MIX Multi-Constituent Urine Control

19480010	Cutoff -50%	5 mL	\$30.00
19480021	Cutoff +50%	5 mL	\$30.00

Advance Diagnostics NZ LTD
 PO BOX 23233
 Hunters Corner
 Auckland
 Tel: 09-2748060
 Fax: 09-2748020
 Email --info@advancediagnostics.co.nz