

- A Handcrafted** in the MDF® artisan tradition by master craftsmen since 1971
- B Clear MDF® ComfortSeal™** eartips provide improved hygiene and comfort
- C Innovative X-configuration** acoustic tubing optimizes acoustic integrity
- D Adult and pediatric** diaphragms and attachments for proper diagnosis
- E Interchangeable** adult, pediatric, and neonatal bell attachments
- F Full-rotation** chestpiece with dual-output acoustic valve stem
- G Extended** non-stick and crack-resistant acoustic tubing
- H Accufit™** adjustable headset for a custom fit
- I Additional** accessories included
- J Latex-free** stethoscope



deluxe sprague rappaport
STETHOSCOPE

MDF® 767X 2-in-1 Deluxe Sprague Rappaport Stethoscope

The MDF® Deluxe Sprague Rappaport Stethoscope with X-configuration dual-lumen acoustic tubing is the most advanced Sprague Rappaport used by healthcare professionals. Similar to the MDF® 767 Sprague Rappaport, it is five stethoscopes in one for a full range of diagnostic applications and is traditionally trusted to detect faint heart sounds and murmurs. The innovative X-configuration tubing eliminates the rubbing of traditional Sprague Rappaport dual tubing that can cause acoustic interference. Healthcare professionals can assess adult, pediatric, and infant patients using this classic, versatile design.