

- A Handcrafted** in the MDF® artisan tradition by master craftsmen since 1971
- B Stress-tested** over 3,000 times and calibrated 5 times before final assembly
- C High-molecular** polymer nylon cuff and case are abrasion and moisture resistant
- D Universal cuff** is compatible with MDF® and other brands of blood pressure systems
- E Certified** 300 mmHg manometer attains accuracy of +/- 3 mmHg without pin stop
- F Chrome-plated** brass screw-type valve delivers precise deflation control rate
- G Accurate** arterial compression with artery indicator label and index range
- H High-density** stick- and crack-resistant inflation bulb and bladder
- I Latex-free** sphygmomanometer



## professional aneroid SPHYGMOMANOMETER

### MDF® 808B Professional Blood Pressure Measurement Device

With its lifetime calibration, the MDF® Professional Aneroid Sphygmomanometer is ideal for medical professionals in all settings. Stress-tested over 3,000 times and calibrated 5 times before final assembly, it is one of the most durable, accurate, and dependable blood pressure measurement devices on the market.