

- A Handcrafted** in the MDF® artisan tradition by master craftsmen since 1971
- B Ergonomic** palm design enables precise single-hand air inflation and release
- C Certified** 300 mmHg manometer attains accuracy of +/- 3 mmHg without pin
- D Universal** cuff is compatible with MDF® and other brands of blood pressure systems
- E High-molecular** polymer nylon cuff and case are abrasion and moisture resistant
- F Chrome-plated** brass screw-type valve delivers precise deflation control rate
- G Stress-tested** over 3,000 times and calibrated 5 times before final assembly
- H Extra-large** high-contrast dial ensures easy and accurate readings
- I High-density** stick- and crack-resistant inflation bulb and bladder
- J Latex-free** sphygmomanometer



medic palm aneroid™
SPHYGMOMANOMETER

MDF® 848XP Professional-Grade Aneroid Sphygmomanometer

The single-hand operated MDF® Medic Palm™ Aneroid Sphygmomanometer features an extra-large faceplate and shock-absorbing gauge guard. With inflation bulb and valve integrated into the gauge, this lightweight super-duty model combines performance, efficiency, and portability. Stress-tested over 3,000 times and calibrated 5 times before final assembly, it is one of the most durable, accurate, and dependable blood pressure measurement devices on the market.