

DRUG OF ABUSE

Cut-off Levels

Class of Drug	Sure Step Test Cup & Cassettes Cut-off level ng/L	Australian Standard AS4308:2008 (Table 1 Immunoassay) ng /L
Morphine (opiates)	300	300
Cannabis metabolites	50	50
Cocaine metabolites	300	300
Benzodiazepines	200	200
Methamphetamines*	300	300

* Please note, the Australian/NZ Standard for urine testing AS4308 shows the 'Class of Drug' as Sympathomimetic amines which includes drug such as Methamphetamines and Amphetamines