

# ACCU-CHEK® Performa

## Product Specifications

Feature	Specification
<b>Test principle</b>	Glucose dehydrogenase with indication by electrochemical alternating current/direct current method
<b>Test range</b>	0.6-33.3 mmol/L
<b>Test period</b>	5 seconds
<b>Ambient temperature during use</b>	6-44°C
<b>Air humidity during use</b>	10-90%
<b>Sample volume</b>	0.6 µL
<b>Haematocrit range</b>	10-70%
<b>Maximum elevation</b>	3,094 meters' above sea level
<b>Test material</b>	Capillary, venous, arterial and neonatal blood, alternative site testing, AST: forearm, upper arm, calf, thigh & thumb
<b>New help and safety functions</b>	<ul style="list-style-type: none"><li>• Under dosage detection</li><li>• Automatic recognition and marking of control solutions</li><li>• Automatic recognition and rejection of expired test strips</li><li>• Date and time preset on delivery</li><li>• Internal back-up battery for maintenance of date and time</li><li>• Long battery life</li></ul>

<b>Calibration</b>	A code key is included in each vial of test strips
<b>Test and hypo alarms</b>	<ul style="list-style-type: none"><li>• Programmable for test reminders</li><li>• Warns when hypoglycemic level is reached</li></ul>

Available at

**Advance Diagnostics NZ Ltd**

PO BOX 217107 Botany Junction Auckland

Auckland

Tel :09-2748060

Fax:09-2748020

[advancediagnostics@gmail.com](mailto:advancediagnostics@gmail.com)

[www.advancediagnostics.co.nz](http://www.advancediagnostics.co.nz)