

# Urine Control With Biochemical Diagnostics

## Detectabuse™ AU/NZ MIX

### For Hand Held Devices



#### Top Reasons to use Detect abuse™ Controls:

- Human urine-based matrix simulates patient samples
- Recommended by the majority of handheld device manufacturers
- 12 drug constituents within one control vial
- Stable for up to 3 years frozen, or refrigerated (unopened)
- Ready to use in 5 mL vial sizes
- Use Detect abuse™ controls to train staff in the accurate use of hand-held devices
- Use Detect abuse™ controls to evaluate the performance of various devices on the market
- Each Detect abuse™ control contains: THC, Benzoylcegonine (COC), PCP, Morphine (OPI), Methamphetamine, Amphetamine, Secobarbital (BAR), Oxazepam (BZO), Methadone, Methaqualone, Propoxyphene, Nortriptyline (TCA).

#### The Detectabuse Advantage:

Detectabuse AU/NZ MIX controls are formulated to work with all known handheld drug-screening devices. We also offer a complete line of controls for immunoassay and GC/MS

#### NEW! AUSTRALIAN-NEW ZEALAND CUTOFF LEVELS

TARGET VALUES	-50%	Cutoff*	+50%
Delta-9-THC-COOH	25	50	75
Benzoylcegonine	150	300	450
Phencyclidine (PCP)	12.5	25	37.5
Morphine	150	300	450
d-Methamphetamine	150	300	450
d-Amphetamine	150	300	450
Secobarbital	150	300	450
Oxazepam	100	200	300
Methadone	150	300	450
Methaqualone	150	300	450
Propoxyphene	150	300	450
Nortriptyline	500	1000	1500

#### Detectabuse AU/NZ MIX Multi-Constituent Urine Control

19480010	Cutoff -50%	5 mL	\$35.00
19480021	Cutoff +50%	5 mL	\$35.00

**Advance Diagnostics NZ Ltd**  
**PO Box 217-107**  
**Botany Junction**  
**Auckland 2164**  
**0800283783( 08002DRUGTEST)**